




SY'AH IMPEX



MAKING PATIENT SAFETY A STANDARD

shifa[®]

PAEDIATRIC INFUSION SET

shifa® PAEDIATRIC INFUSION SET


Time management is a key challenging factor for health care professionals to perform variety of tasks effectively and efficiently for patient's safety and satisfaction.

Shifa I.V. Burette contains a metered small volume chamber designed to limit the amount of solution available to the patient, usually for safety, but also an intermittent infusion of medication that can be mixed immediately before infusion. Primarily required for pediatric and neonatal patients.

Suitable for gravity infusions, the use of the burette allows fine control of fluid volume and avoid inadvertent over transfusion.

Benefits at a glance

- Blister packed ensure sterilization.
- Ethylene Oxide (EO) sterilized.
- Delivery of drops at exact measurement.
- Gravity feed -no external power requirements.
- Functions with minimal monitoring.
- Simple to use, easy staff training.
- Soft material chamber can accept push.
- Designed to meet international burette standard ISO 8536-5:2004.
- CE-marked and manufactured in accordance with ISO 13485:2003 certified facility.



A product Of SY'AH IMPEX
1-6/15, Sector 5, Korangi Industrial Area, Karachi.
Ph: (92 21) 35122997-98, Fax: (92 21) 35065455
email: info@syahimpex.com
f shifabrands | t shifabrands



SY'AH IMPEX provides this brochure for informational purpose only. It is not intended to be a substitute for professional medical advice, diagnosis or treatment.