



SY'AH IMPEX



SETTING THE STANDARD IN PROTECTION

**shifa**<sup>®</sup>

SURGICAL GLOVES

# shifa<sup>®</sup> SURGICAL GLOVES

SHIFA surgical gloves act as a barrier to bacteria and infectious viral pathogens widely circulating in surgical examination room and provide safety to HCP during surgical procedures.

## LIGHTLY POWDERED

Allows easy slipping on your hand.

## BARRIER TO ALLERGIES

Low protein content that reduces the risk of allergy.

## HIGHER TENSILE STRENGTH

high flexibility guarantees high levels of elasticity, minimizes bursting when putting on.

## FULLY SECURE

Rolling collar for easy putting on the prevention of back-roll.

## ANATOMICAL SHAPE

Fully anatomical shape gives a natural feeling of comfort without wrinkling or creasing and reduces hand-fatigue during longer period of use.

## MICRO- ROUGHENED SURFACE

Micro-roughened surface that allows easier gripping and handling of instruments in both damp and dry conditions.

## Functional Benefits

- High strength and elasticity of natural high grade latex.
- Intrinsic softness provides superior comfort and natural fit.
- Extra length provides addition protection from surgical debris
- Micro roughened surface gives superiors comfort and natural fit
- Fully anatomical design to reduce hand fatigue during prolonged surgical procedures
- Custom- moulded in various size to fit end user's hand with added comfort and improved fit

## Clinical specifications

- Reducing potential allergic reactions.
- Moisture contents is below 0.8% per glove.
- Sterilized by gamma radiation /ETO process.
- Shelf life is 5 years from date of manufacturing.
- Colour, natural to white and uniform allowing for gradual shift due to thickness variation bioburden and sterility tested.



Dimension	Standards		
	Size	Shifa Gloves	ASTM D3577
Width (mm)	5.5	71 +/-5	(06)
	6.0	77 +/-5	70 +/-6
	6.5	83 +/-5	76 +/-6
	7.0	89 +/-5	83 +/-6
	7.5	95 +/-5	89 +/-6
	8.0	102 +/-5	95 +/-6
	8.5	108 +/-5	102 +/-6
9.0	115 +/-5	108 +/-6	
Lenght(mm)	For all sizes	min 280	114 +/-6
Thickness (mm) Single wall			
Finger	For all sizes	min 0.17	min 0.10
Plam	For all sizes	min 0.12	min 0.10

Dimension	Standards	
	Shifa Gloves	ASTM D3577 (06)
Elongation at break (%)		
Befor aging	min 750	min 750
After aging	min 560	min 560
Tensile Streght (MPa)		
Befor aging	min 24.0	min 24.0
After aging	min 18.0	min 18.0



A product Of SY'AH IMPEX  
1-6/15, Sector 5, Korangi Industrial Area, Karachi.  
Ph: (92 21) 35122997-98, Fax: (92 21) 35065455  
email: info@syahimpex.com  
Made in Malaysia

f shifabrands | t shifabrands



SY'AH IMPEX provides this brochure for informational purpose only. It is not intended to be a substitute for professional medical advice, diagnosis or treatment.